

BED SETTEE WITH 2 ARMS MULTY

Designer

Claude Brisson

Technical characteristics

Multy offers an intermediate position halfway between settee and bed: a chaise longue with ample leg support. Frame in multi-ply and panels of particleboard, with steel mechanism and birch slats (thickness 9 mm). Seat in either polyether foam (28kg/m3 - 3.0kPa) or 9 cm thick high resilience polyurethane Bultex foam, density 38 kg/ m3. Back in polyether foam (28kg/m³ - 3.0kPa).Back in polyether foam (28kg/m³ - 3.0kPa). Lumbar cushions are supplied as follows: 2 x W60 for small bed-settee without arms, and 2 x W78 for large bed-settee without arms; 1 x W78 for small bedsettee with arms, and 2 x W60 for large bed-settee with arms. Please take these details into consideration when ordering additional/subsequent lumbars. The armrests (polyether foam 28 kg/m³ - 3.0 kPa + polyester quilting 200 g/m²) or bedside tables (MDF top diameter 43 cm - thickness 2 cm) pivot to allow the bedsettee to open. MULTY is easy to operate, and offers excellent comfort both as a bed (with polyether mattress size: 153x194x8, 119x194x8) and as a settee. When unfolded, its depth is 200 cm for the bed-settee W 155. Mattress has an ecru-coloured cover in 50% polyester/50% cotton (dry cleanable), which has been treated with "HealthGuard" (antibacterial &



DIMENSIONS

fungicidal). A foam sheet is supplied to place between the slats and the seat when folded. For instructions on use, please see the card which comes with the bed. MULTY has been specially designed for use as a bed on an everyday basis. Optional extras: lumbar cushions and arms. Covers of lumbar cushions and armrests are removable.

More information on

www.cinna.fr © Ligne Roset 2024