

BED SETTEE MULTY

Designer

Claude Brisson



Dimensions

Technical characteristics

Frame in multi-ply and panels of particleboard, with steel mechanism and birch slats (thickness 9 mm). Seat in either polyether foam (28kg/m³ - 3.0kPa) or 9 cm thick high resilience polyurethane Bultex foam, density 38 kg/ m³. Back in polyether foam (28kg/m³ - 3.0kPa). Back in polyether foam (28kg/m³ - 3.0kPa). Lumbar cushions are supplied as follows: 2 x W60 for small bed-settee without arms, and 2 x W78 for large bed-settee without arms; 1 x W78 for small bed-settee with arms, and 2 x W60 for large bed-settee with arms. Please take these details into consideration when ordering additional/subsequent lumbar. The armrests (polyether foam 28 kg/m³ - 3.0 kPa + polyester quilting 200 g/m²) or bedside tables (MDF top diameter 43 cm - thickness 2 cm) pivot to allow the bed-settee to open. MULTY is easy to operate, and offers excellent comfort both as a bed (with polyether mattress size : 153x194x8, 119x194x8) and as a settee. When unfolded, its depth is 200 cm for the bed-settee W 155. Mattress has an ecru-coloured cover in 50% polyester/50% cotton (dry cleanable), which has been treated with "HealthGuard" (antibacterial & fungicidal). A foam sheet is supplied to place between the slats and the seat when folded. For instructions on use, please see the card which comes with the bed. MULTY has been specially designed for use as a bed on an everyday basis. Optional extras: lumbar cushions and arms. Covers of lumbar

cushions and armrests are removable. Ready to sleep bed linen includes 2 sheets (top and bottom) that are attached at the foot by a seam. The primary intention is as practical linen for sofa-beds. WASHABLE AT 60deg.C.